

# TREATMENT GUIDE

## Richie Brace® Products, Indications and Casting



### Richie Brace - Standard

- Mild to Moderate PTTD
- Chronic Ankle Instability

Casting:  
- Ankle casting sock



### Medial Arch Suspender

- Moderate to Severe PTTD
- Stage II or III PTTD with subluxed T-N Joint

Casting:  
- Ankle casting sock



### Richie Solid AFO

- Severe Dropfoot - Spasticity
- Charcot Arthropathy
- Dropfoot with unstable knee

Casting:  
- Full Leg (Bermuda) Casting Sock



### Lateral Arch Suspender

- Peroneal tendinopathy fixed varus deformity of hindfoot/ankle
- Severe Lateral ankle Instability

Casting:  
- Ankle casting sock



### Richie OTC Dynamic Assist

- Paralysis or Weakness of Ankle Muscles
- Peroneal or Tibialis Anterior Tendinopathy

Casting:  
- No Casting Necessary  
Order by Left/Right  
S, M, L and XL



### Richie Brace - Restricted Hinge

- DJD of Ankle or rearfoot
- Dropfoot - Mild
- Dropfoot with equinus
- Dropfoot with spasticity
- Stable knee

Casting:  
- Ankle casting sock



### Richie Gauntlet

- Rigid, non-reducible Adult Acquire Flatfoot (Stage III & IV)
- Severe DJD or Deformity of hindfoot
- Charcot Arthropathy

Casting:  
- Mid Leg Casting Sock



### Richie OTC Ankle Brace

- Prefabricated Ankle Brace for:
- Acute Ankle sprain
- Tendinitis of Ankle, Interim Brace
- Before Custom AFO treatment

Casting:  
- No Casting Necessary  
Order by Left/Right  
S, M, L and XL



### Dynamic Assist

- Dropfoot without Equinus
- No spasticity
- Stable knee (must have all 3 above)
- Post CVA
- Peroneal nerve injury
- Post Polio

Casting:  
- Ankle casting sock



### California AFO

- Gauntlet brace without laces
- Severe Deformity: Stage IV PTTD
- Charcot Deformity
- Severe DJD of Ankle or Hindfoot

Casting:  
- Mid Leg Casting Sock