

Care Instructions

1. The Richie Brace® is a custom ankle foot orthosis which has several component parts, including:

- A plastic foot plate and two plastic leg uprights
- Soft interface padding along the leg uprights
- Foam top cover on the foot plate
- Velcro strap closures

2. Several other component parts may be incorporated into your brace including:

- Arch suspender strap and cork arch flap
- Ankle hinge joints made of metal or plastic
- Additional padding under the toes and ball of the foot
- Leather boot closure (California AFO)

3. The following care instructions pertain to these component parts:

- Never immerse the brace into water
- Never expose the brace to heat
- Do not attempt to apply glue or additional padding to any part of the brace
- Do not apply any cleaning agents to any portion of the brace

4. Almost all component parts of the Richie Brace® can be replaced or refurbished after extended periods of use. The padding, Velcro straps and hinges can be replaced at the discretion of your brace provider.

